



DEVICE DETAILS

NAME OF DEVICE	Vircast 3D Comfort Fit Medical Face Mask
ESTABLISHMENT NAME	VIRCAST MEDIC SDN BHD
ROLE OF ESTABLISHMENT	MANUFACTURER
REGISTRATION NO	GA2984124-175552
BRAND NAME	Vircast Medic
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	The breathable waterproof mask made of non-woven fabric (outer) and (inner) melt blown fabric filtering layer that effectively waterproof and isolate droplets. This mask entirely made up of multiple filter material layers that can give the optimum protection and comfort.
DEVICE INTENDED PURPOSE	This face mask helps to block dusts, germs, dirt, and pollen while also absorbing excess moisture. It reduces the number of germs that the wearer releases while protecting from being infected.
VALIDITY DATE OF REGISTRATION	03/07/2024 - 02/07/2029

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Vircast 3D Comfort Fit Medical Face Mask	VM3D002
2	Vircast 3D Comfort Fit Medical Face Mask	VM3D003
3	Vircast 3D Comfort Fit Medical Face Mask	VM3D004
4	Vircast 3D Comfort Fit Medical Face Mask	VM3D005
5	Vircast 3D Comfort Fit Medical Face Mask 0-3	VM3D006
6	Vircast 3D Comfort Fit Medical Face Mask 3-6	VM3D007
7	Vircast 3D Comfort Fit Medical Face Mask	VM3DHJ001