



DEVICE DETAILS

NAME OF DEVICE	FITSCAN 365 PULSE OXIMETER
ESTABLISHMENT NAME	MAXWELL MEDICAL SUPPLIES SDN BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GB1845921-81484
BRAND NAME	FITSCAN 365
MEDICAL DEVICE CATEGORY	MD 1302 - Monitoring devices of vital physiological parameters
DEVICE GROUPING TYPE	SINGLE
DEVICE DESCRIPTION	FitScan 365 Pulse Oximeter is used to measure arterial oxygen saturation (%SpO ₂) of hemoglobin and pulse rate, an important indicator of respiratory function, using a pair of sensors placed on a thin part of the patient's fingertip.
DEVICE INTENDED PURPOSE	FitScan 365 Pulse Oximeter is intended for spot checking of functional oxygen saturation of arterial hemoglobin (SpO ₂) and pulse rate. This portable device is indicated for use in adult patients in clinical institution and home environments.
VALIDITY DATE OF REGISTRATION	21/12/2021 - 20/12/2026

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	FITSCAN 365 PULSE OXIMETER	0785104813816