



DEVICE DETAILS

NAME OF DEVICE	FITSCAN 365 PULSE OXIMETER
ESTABLISHMENT NAME	MAXWELL MEDICAL SUPPLIES SDN BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GB1845921-81484
BRAND NAME	FITSCAN 365
MEDICAL DEVICE CATEGORY	MD 1302 - Monitoring devices of vital physiological parameters
DEVICE GROUPING TYPE	SINGLE
DEVICE DESCRIPTION	FitScan 365 Pulse Oximeter is used to measure arterial oxygen saturation (%SpO2) of hemoglobin and pulse rate, an important indicator of respiratory function, using a pair of sensors placed on a thin part of the patient's fingertip.
DEVICE INTENDED PURPOSE	FitScan 365 Pulse Oximeter is intended for spot checking of functional oxygen saturation of arterial hemoglobin (SpO2) and pulse rate. This portable device is indicated for use in adult patients in clinical institution and home environments.
VALIDITY DATE OF REGISTRATION	21/12/2021 - 20/12/2026

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	FITSCAN 365 PULSE OXIMETER	0785104813816