




## DEVICE DETAILS

NAME OF DEVICE	MEDICOS M99 LITE FLIP SURGICAL FACE MASK
ESTABLISHMENT NAME	CREATIVE CONTRACT (M) SDN BHD
ROLE OF ESTABLISHMENT	MANUFACTURER
REGISTRATION NO	GA7828624-176362
BRAND NAME	MEDICOS
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	- 5-ply Super Soft PP Non-Woven - Bacterial Filtration Efficiency $\geq 99\%$ - Particle Filtration Efficiency $0.1\mu \geq 99\%$ - Fluid Resistant
DEVICE INTENDED PURPOSE	Used inside the operating room or within other sterile procedure areas to protect the patient environment from contamination. It is also intended to protect the wearer against bacteria, viruses, particles and splashes from potentially contaminated bio-liquids.
VALIDITY DATE OF REGISTRATION	12/07/2024 - 11/07/2029

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	M99 LITE FLIP SURGICAL FACE MASK 	9800
2	M99 LITE FLIP SURGICAL FACE MASK	9801
3	M99 LITE FLIP SURGICAL FACE MASK	9802
4	M99 LITE FLIP SURGICAL FACE MASK	9803
5	M99 LITE FLIP SURGICAL FACE MASK	9804
6	M99 LITE FLIP SURGICAL FACE MASK	9805
7	M99 LITE FLIP SURGICAL FACE MASK	9806
8	M99 LITE FLIP SURGICAL FACE MASK	9807
9	M99 LITE FLIP SURGICAL FACE MASK	9808
10	M99 LITE FLIP SURGICAL FACE MASK	9809