



DEVICE DETAILS

NAME OF DEVICE	CHILINE PULSE OXIMETER MONITORING SYSTEM (WITH OR WITHOUT MONITORING RESPIRATION)
ESTABLISHMENT NAME	WELLSPRING SOLUTIONS SDN. BHD.
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GC6331224-177056
BRAND NAME	CHILINE/ S68 AND S68 LITE
MEDICAL DEVICE CATEGORY	MD 1100 - GENERAL ACTIVE MEDICAL DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	This product is a monitoring system for blood oxygen level and respiration. It can simultaneously measure blood oxygen level (SpO ₂), heart rate and respiratory rate, and can record and display results.
DEVICE INTENDED PURPOSE	This product uses optical sensor, a non-invasive technology that is suitable for adults, children and infants. Upon recommendations from their physicians, users can operate this product on their own to measure single or continuous blood oxygen level, heart rate and respiratory rate. The results will be displayed in a mobile App, serving as a reference for physicians for diagnosis. Note: Model S68 Lite does not have a respiratory monitoring feature.
VALIDITY DATE OF REGISTRATION	20/07/2024 - 19/07/2029

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Chiline Pulse Oximeter Monitoring System (with or without Monitoring Respiration)	S68
2	Chiline Pulse Oximeter Monitoring System (with or without Monitoring Respiration)	S68 Lite
3	Storage case	N/A
4	Charging cable	N/A
5	Protective cover	N/A
6	Instruction manual	N/A
7	L-Size SpO ₂ Probe	PA68
8	M-Size SpO ₂ Probe	PM68
9	S-Size SpO ₂ Probe	PC68

NO	NAME OF DEVICE	IDENTIFIER
10	XS-Size SpO2 Probe	PB68
11	Ring SpO2 Probe	PR68
12	Chiline HomeCare App	Version: V3



MEDICAL DEVICE AUTHORITY