




DEVICE DETAILS

NAME OF DEVICE	Medical Face Mask
ESTABLISHMENT NAME	H&A PRODUCTS SDN BHD
ROLE OF ESTABLISHMENT	MANUFACTURER
REGISTRATION NO	GA2664324-177457
BRAND NAME	SHIRUDO Medical Face Mask
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	SINGLE
DEVICE DESCRIPTION	SHIRUDO 3-Ply Medical Face Mask is made of 3 layers of high quality non-woven fabric , meeting the requirements of ASTM F2100 and BS EN 14683. It achieves Bacterial Filtration Efficiency (BFE) of more than 99% and Sub-micron Particulate Filtration Efficiency (PFE) of more than 98%
DEVICE INTENDED PURPOSE	SHIRUDO 3-Ply Medical Face Mask is designed for single use and it provides optimum protection against transfer of microorganisms, blood fluid and airborne particulate
VALIDITY DATE OF REGISTRATION	30/07/2024 - 29/07/2029

LIST OF DEVICE

NO	NAME OF DEVICE 	IDENTIFIER
1	Medical Face Mask	FMFF3L02AS