




DEVICE DETAILS

NAME OF DEVICE	SAFETYWARE 3-Ply Medical Face Mask Ear Loop & Head Loop
ESTABLISHMENT NAME	SAFETYWARE SDN. BHD.
ROLE OF ESTABLISHMENT	MANUFACTURER
REGISTRATION NO	GA9771324-180795
BRAND NAME	SAFETYWARE
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	A 3-ply face mask is a medical device designed for respiratory protection. Its key characteristics include three layers of material – typically a non-woven outer layer, a melt-blown middle layer for filtration, and a soft inner layer for comfort. The materials used are often polypropylene. The design aims to provide effective filtration of airborne particles, offering a barrier against respiratory droplets.
DEVICE INTENDED PURPOSE	Provide a physical barrier that helps prevent the spread of respiratory droplets.
VALIDITY DATE OF REGISTRATION	26/08/2024 - 25/08/2029

LIST OF DEVICE

NO	NAME OF DEVICE 	IDENTIFIER
1	SAFETYWARE 3-Ply Medical Face Mask (Earloop)	1302EL-BLU
2	SAFETYWARE 3-Ply Medical Face Mask (Headloop)	1302HL-BLU