



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	SAFETYWARE 3-Ply Medical Face Mask Ear Loop & Head Loop
<b>ESTABLISHMENT NAME</b>	SAFETYWARE SDN. BHD.
<b>ROLE OF ESTABLISHMENT</b>	MANUFACTURER
<b>REGISTRATION NO</b>	GA9771324-180795
<b>BRAND NAME</b>	SAFETYWARE
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	FAMILY
<b>DEVICE DESCRIPTION</b>	A 3-ply face mask is a medical device designed for respiratory protection. Its key characteristics include three layers of material – typically a non-woven outer layer, a melt-blown middle layer for filtration, and a soft inner layer for comfort. The materials used are often polypropylene. The design aims to provide effective filtration of airborne particles, offering a barrier against respiratory droplets.
<b>DEVICE INTENDED PURPOSE</b>	Provide a physical barrier that helps prevent the spread of respiratory droplets.
<b>VALIDITY DATE OF REGISTRATION</b>	26/08/2024 - 25/08/2029

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	SAFETYWARE 3-Ply Medical Face Mask (Earloop)	1302EL-BLU
2	SAFETYWARE 3-Ply Medical Face Mask (Headloop)	1302HL-BLU