



DEVICE DETAILS

NAME OF DEVICE	KF94 Medical Mask
ESTABLISHMENT NAME	VIRCAST MEDIC SDN BHD
ROLE OF ESTABLISHMENT	MANUFACTURER
REGISTRATION NO	GA10661524-187952
BRAND NAME	Vircast Medic
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	The medical face mask generally composed of a filter layer that is placed, bonded or molded between layers of fabric. The disposal medical face mask is mainly composed of three layers of non-woven fabric with a pair of elastic ear loop and a piece of nose piece in which being assembled using ultrasonic sewing technology.
DEVICE INTENDED PURPOSE	The main intended use of medical face masks is to protect the patient from infective agents and help reduce exposure of your saliva n respiratory secretions to others. Medical face masks may also be intended to be worn by patients and healthcare professional to reduce the risk of spread of infections, particularly in epidemic or pandemic situations.
VALIDITY DATE OF REGISTRATION	01/11/2024 - 31/10/2029

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	KF94 MEDICAL FACE MASK	KF0001
2	KF94 MEDICAL FACE MASK	KF0002
3	KF94 MEDICAL FACE MASK - CANDY	KF0003
4	KF94 MEDICAL FACE MASK - GLAMOROUS	KF0004
5	KF94 KIDS MEDICAL FACE MASK - MIX	KF0005
6	KF94 KIDS MEDICAL FACE MASK	KF0006
7	KF94 MEDICAL FACE MASK - SHARMEELA	KF0007
8	KF94 MEDICAL FACE MASK - HIJABIS	KF0008
9	KF94 MEDICAL FACE MASK - MIX	KF0009

NO	NAME OF DEVICE	IDENTIFIER
10	KF94 MEDICAL FACE MASK	KF0010



MEDICAL DEVICE AUTHORITY