

DEVICE DETAILS

| | |
|-------------------------------|---|
| NAME OF DEVICE | KF94 Medical Mask |
| ESTABLISHMENT NAME | VIRCAST MEDIC SDN BHD |
| ROLE OF ESTABLISHMENT | MANUFACTURER |
| REGISTRATION NO | GA10661524-187952 |
| BRAND NAME | Vircast Medic |
| MEDICAL DEVICE CATEGORY | SINGLE-USE DEVICES |
| DEVICE GROUPING TYPE | FAMILY |
| DEVICE DESCRIPTION | The medical face mask generally composed of a filter layer that is placed, bonded or molded between layers of fabric. The disposal medical face mask is mainly composed of three layers of non-woven fabric with a pair of elastic ear loop and a piece of nose piece in which being assembled using ultrasonic sewing technology. |
| DEVICE INTENDED PURPOSE | The main intended use of medical face masks is to protect the patient from infective agents and help reduce exposure of your saliva n respiratory secretions to others. Medical face masks may also be intended to be worn by patients and healthcare professional to reduce the risk of spread of infections, particularly in epidemic or pandemic situations. |
| VALIDITY DATE OF REGISTRATION | 01/11/2024 - 31/10/2029 |

LIST OF DEVICE

| NO | NAME OF DEVICE | IDENTIFIER |
|----|------------------------------------|------------|
| 1 | KF94 MEDICAL FACE MASK | KF0001 |
| 2 | KF94 MEDICAL FACE MASK | KF0002 |
| 3 | KF94 MEDICAL FACE MASK - CANDY | KF0003 |
| 4 | KF94 MEDICAL FACE MASK - GLAMOROUS | KF0004 |
| 5 | KF94 KIDS MEDICAL FACE MASK - MIX | KF0005 |
| 6 | KF94 KIDS MEDICAL FACE MASK | KF0006 |
| 7 | KF94 MEDICAL FACE MASK - SHARMEELA | KF0007 |
| 8 | KF94 MEDICAL FACE MASK - HIJABIS | KF0008 |
| 9 | KF94 MEDICAL FACE MASK - MIX | KF0009 |

| NO | NAME OF DEVICE | IDENTIFIER |
|----|------------------------|------------|
| 10 | KF94 MEDICAL FACE MASK | KF0010 |



MEDICAL DEVICE AUTHORITY