

DEVICE DETAILS		
NAME OF DEVICE	JUMPER PULSE OXIMETER JPD-500D, JPD-500A JPD-500C JPD-500E JPD-500F JPD-500B JPD-500G JPD-500H	
ESTABLISHMENT NAME	META PHARMA SDN BHD	
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE	
REGISTRATION NO	GC10853124-162786	
BRAND NAME	JUMPER	
MEDICAL DEVICE CATEGORY	MD 0104 - Non-active medical devices with measuring function	
DEVICE GROUPING TYPE	FAMILY	
DEVICE DESCRIPTION	The fingertip pulse Oximeter features small size, low power consumption, convenient operation and portability. It is only necessary for a user to put one of his fingers into the fingertip photoelectric sensor for diagnosis, and a display screen will show oxygen saturation.	
DEVICE INTENDED PURPOSE	The Fingertip Pulse Oximeter is non-invasive device intended for spot- checking of functional oxygen saturation of arterial hemoglobin (Spo2) and pulse rate. The portable fingertip device is indicated for adult and pediatric patients in home and hospital environments.	
VALIDITY DATE OF REGISTRATION	06/02/2024 - 05/02/2029	

LIST OF DEVICE		
NO	NAME OF DEVICE	IDENTIFIER
1	JPD-500D	JPD-500D