



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	Electric Wheelchair Lightweight (Airline Approved)
<b>ESTABLISHMENT NAME</b>	FRESCO COCOA SUPPLY PLT
<b>ROLE OF ESTABLISHMENT</b>	MANUFACTURER
<b>REGISTRATION NO</b>	GA6287724-188892
<b>BRAND NAME</b>	FRESCO bike
<b>MEDICAL DEVICE CATEGORY</b>	ASSISTIVE PRODUCTS FOR PERSONS WITH DISABILITY
<b>DEVICE GROUPING TYPE</b>	SINGLE
<b>DEVICE DESCRIPTION</b>	this electric wheelchair come with push handle, joystick, armrest, armrest lifting switch, 12" rear (honeycomb) wheel, battery, 8" front (honeycomb), anti reverse wheel, manual/power switch, and locking path
<b>DEVICE INTENDED PURPOSE</b>	USED BY PEOPLE FOR WHOM WALKING DIFFICULT OR IMPOSSIBLE DUE TO ILLNESS,INJURY OR DISABILITY.
<b>VALIDITY DATE OF REGISTRATION</b>	08/11/2024 - 07/11/2029

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Electric Wheelchair Lightweight (Airline Approved)	73001