




## DEVICE DETAILS

NAME OF DEVICE	Disposable Face Mask
ESTABLISHMENT NAME	ADVENTA HEALTHCARE SDN. BHD.
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GA2281022-91384
BRAND NAME	Connecx Respinecx
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	This disposable face mask is made of spun polyethylene + melt-blown polyethylene + spun polyethylene.
DEVICE INTENDED PURPOSE	Disposable Face Mask is used to cover the mouth and nose to serves as a barrier for contaminants from getting released from the mouth of a sick person and spread to the environment or to be inhaled by people in the public.
VALIDITY DATE OF REGISTRATION	27/04/2022 - 26/04/2027

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Disposable Face Mask 	UM120-B
2	Disposable Face Mask	UM120-G
3	Disposable Face Mask	UM120P-B
4	Disposable Face Mask	UM120P-G
5	Disposable Face Mask	UM121-B
6	Disposable Face Mask	UM121-G
7	Disposable Face Mask	UM511-B
8	Disposable Face Mask	UM511-G