




DEVICE DETAILS

NAME OF DEVICE	Surgical Face Mask
ESTABLISHMENT NAME	QUANTUM UPSTREAM SDN.BHD.
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GA6244825-204592
BRAND NAME	HospiGuard
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	SINGLE
DEVICE DESCRIPTION	Surgical Face Mask (Ear Loop) is made from 2 layers of non-woven fabric and 1 layer of melt-blown filter fabric with high BFE over 95%, hypoallergenic and very low resistance to breath
DEVICE INTENDED PURPOSE	Surgical Face Mask is intended for used as a barrier to prevent air-borne transfer of bacteria and particles. This is intended for used in operating theater.
VALIDITY DATE OF REGISTRATION	09/06/2025 - 08/06/2030

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Surgical Face Mask 	MF050101