



DEVICE DETAILS

NAME OF DEVICE	Anti-Choking Device
ESTABLISHMENT NAME	MEDIHOMEAID PLT
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GA3263925-214726
BRAND NAME	Sonmol
MEDICAL DEVICE CATEGORY	REUSABLE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	The device is a portable anti-choking device designed to assist in choking by using negative pressure. Its working principle involves a spring-loaded piston mechanism: pressing the top piston compresses the internal spring and stores energy. When the release button is pressed, the spring expands, driving the piston upward and creating a negative pressure chamber between the piston and the outer shell. When the mask is properly fitted over the user's mouth and nose, this negative pressure helps extract airway obstructions safely and effectively.
DEVICE INTENDED PURPOSE	It is a non-invasive, reusable device designed to assist in choking emergency.
VALIDITY DATE OF REGISTRATION	24/10/2025 - 23/10/2030

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Sonmol Anti-Choking Device	MFAC-08VS
2	Sonmol Anti-Choking Device	MFAC-08VC