



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	BLOOD PRESSURE MONITOR MEDISENSE X9
<b>ESTABLISHMENT NAME</b>	PFLEGE MEDICAL SUPPLIES SDN. BHD.
<b>ROLE OF ESTABLISHMENT</b>	MANUFACTURER
<b>REGISTRATION NO</b>	GB10255725-216438
<b>BRAND NAME</b>	PROFIX
<b>MEDICAL DEVICE CATEGORY</b>	MD 1100 - GENERAL ACTIVE MEDICAL DEVICES
<b>DEVICE GROUPING TYPE</b>	SINGLE
<b>DEVICE DESCRIPTION</b>	PROFIX Blood Pressure Monitor Medisense X9 consists of eight main components and accessories, including LCD, Keys, Pump, Cuff, MCU, Transducer, Energy Source and Case. This is an automatic non-invasive blood pressure monitor which can be driven by 4 AA batteries or type-C USB (optional). It uses an inflatable cuff which is wrapped around the patient's upper arm to measure the systolic and diastolic blood pressure as well as the pulse rate of adult at household, not for neonate or pregnancy.
<b>DEVICE INTENDED PURPOSE</b>	PROFIX Blood Pressure Monitor Medisense X9 is intended purposed to measure the systolic and diastolic blood pressure and pulse rate of an adult individual by using a non-invasive technique in which an inflatable cuff is wrapped around the upper arm.
<b>VALIDITY DATE OF REGISTRATION</b>	13/11/2025 - 12/11/2030

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	BLOOD PRESSURE MONITOR MEDISENSE X9	MEDISENSE X9