



DEVICE DETAILS

NAME OF DEVICE	GM85 (C-TYPE, F-TYPE, FIT-TYPE)
ESTABLISHMENT NAME	LAC MEDICAL SUPPLIES SDN BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GC5394423-134640
BRAND NAME	SAMSUNG ELECTRONIC CO.LTD.
MEDICAL DEVICE CATEGORY	MD 1201 - Imaging devices utilising ionizing radiation
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	The GM85 digital X-ray imaging system consists of the built-in high voltage generator, X-ray tube, control panel, mobile system for movement, detector and the S-station (operation software).
DEVICE INTENDED PURPOSE	The GM85 Digital X-ray Imaging System is intended for use in generating radiographic images of human anatomy by a qualified/ trained doctor or technician. This device is not intended for mammographic applications.
VALIDITY DATE OF REGISTRATION	08/06/2023 - 07/06/2028

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
No results found.		