



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	LOMOCHA N95 Medical Protective Face Mask
<b>ESTABLISHMENT NAME</b>	BLUE ORIGIN SUPPLY CHAIN SDN. BHD.
<b>ROLE OF ESTABLISHMENT</b>	AUTHORIZED REPRESENTATIVE
<b>REGISTRATION NO</b>	GA4741826-222164
<b>BRAND NAME</b>	N95 Medical Protective Face Mask
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	SINGLE
<b>DEVICE DESCRIPTION</b>	The product is composed of mask body, nose clip and mask band. The mask body is made of non-woven fabric and melt blown fabric. The nose clip is made of flexible materials. The mask band is made of elastic fabric. The product can be worn by all kinds of clinical personnel during non-invasive operation to cover the mouth, nose and lower jaw of the user so as to provide a certain physical barrier to prevent the direct penetration of pathogenic microorganisms and particles.
<b>DEVICE INTENDED PURPOSE</b>	The product can be worn by all kinds of clinical personnel during non-invasive operation to cover the mouth, nose and lower jaw of the user so as to provide a certain physical barrier to prevent the direct penetration of pathogenic microorganisms and particles.
<b>VALIDITY DATE OF REGISTRATION</b>	14/01/2026 - 13/01/2031

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
No results found.		