



DEVICE DETAILS

NAME OF DEVICE	ULTRAFORMER III
ESTABLISHMENT NAME	ELOGIO ASIA PACIFIC SDN BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GB7557426-226700
BRAND NAME	CLASSYS
MEDICAL DEVICE CATEGORY	MD 1100 - GENERAL ACTIVE MEDICAL DEVICES
DEVICE GROUPING TYPE	SYSTEM
DEVICE DESCRIPTION	The ULTRAFORMER III is a non-invasive equipment. It is typically computer-controlled and is capable of producing localized heat within tissues or organs. Energy delivered to the patient is via an externally-mounted transducer. It is a safe and effective device for the treatment of dermatologic diseases.
DEVICE INTENDED PURPOSE	A non-invasive therapeutic device for the treatment of ptosis, axillary hyperhidrosis and correction of skin laxity using the principle of coagulation and necrosis.
VALIDITY DATE OF REGISTRATION	25/02/2026 - 24/02/2031

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	ULTRAFORMER III (Main Body)	3217-00025-601
2	Hand-piece ASSY	UFH000002
3	L4-4.5	UFC344501
4	L7-3.0	UFC373001
5	L7-1.5	UFC371501
6	MF1	UFC3MF103
7	MF2	UFC3MF203
8	MF3	UFC3MF304
9	MF5	UFC3MF504
10	MF6	UFC3MF604
11	MF9	UFC3MF904

NO	NAME OF DEVICE	IDENTIFIER
12	MF13	UFC3MFD04
13	LCD Assembly (Main Body)	K555-10402-001
14	Tocuh Panel	K55-10402-001
15	LCD Module (Hand-piece)	K550-0220H-001
16	Fuse	G650-25010-001
17	Key switch	G010-00001-001
18	AC Power Cord	5620-00001-001
19	AC Inlet	7438-00170-001
20	Relay	D500-00001-001
21	System SMPS	G210-06312-003
22	12V SMPS	G210-00512-001
23	RF SMPS	G210-10848-003



MEDICAL DEVICE AUTHORITY