



DEVICE DETAILS

NAME OF DEVICE	PULSE OXIMETER
ESTABLISHMENT NAME	MEDISANA MALAYSIA SDN BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GC5975226-232412
BRAND NAME	MEDISANA
MEDICAL DEVICE CATEGORY	MD 1100 - GENERAL ACTIVE MEDICAL DEVICES
DEVICE GROUPING TYPE	SINGLE
DEVICE DESCRIPTION	Medisana Pulse Oximeter (PM 180) determines SpO2 by measuring the absorption of red & infrared light passing through perfused tissue. Changes in absorption caused by the pulsation of blood in the vascular bed are used to determine SpO2 reading and pulse rate.
DEVICE INTENDED PURPOSE	The Fingertip Pulse Oximeter are intended for measuring functional oxygen saturation of arterial hemoglobin (SpO2) and pulse rate for adults as non-invasive spot checking in home environment. It is designed for fingers between 0.8cm and 2.5cm (0.3 inches to 1 inches) and for patients during no-motion condition.
VALIDITY DATE OF REGISTRATION	13/05/2026 - 12/05/2031

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
No results found.		