

DEVICE DETAILS

NAME OF DEVICE	TIGHTAN PRO
ESTABLISHMENT NAME	AESTHETIC LEADER SDN BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GB5232326-233211
BRAND NAME	TIGHTAN PRO
MEDICAL DEVICE CATEGORY	MD 1100 - GENERAL ACTIVE MEDICAL DEVICES
DEVICE GROUPING TYPE	SYSTEM
DEVICE DESCRIPTION	The Tightan Pro is a focused ultrasound stimulator system. The ultrasonic frequency of the device is 2.0 MHz, 4.0 MHz, 5.5 MHz and 7.0 MHz and the energy is up to 3 J. The Tightan Pro consists of a main unit, handpieces, cartridges and a power cable.
DEVICE INTENDED PURPOSE	The Tightan Pro is intended for use of non-invasive coagulation of soft tissue.
VALIDITY DATE OF REGISTRATION	04/06/2026 - 03/06/2031

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
No results found.		