

DEVICE DETAILS

NAME OF DEVICE	BUZUD PULSE OXIMETER
ESTABLISHMENT NAME	BUZUD MALAYSIA SDN. BHD.
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GB6901926-239654
BRAND NAME	BUZUD
MEDICAL DEVICE CATEGORY	MD 1100 - GENERAL ACTIVE MEDICAL DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	SpO2 stands for peripheral capillary oxygen saturation. Oxygen saturation is defined as the ratio of oxyhemoglobin (HbO2) to the total concentration of hemoglobin (i.e. Oxyhemoglobin + reduced hemoglobin) present in the blood. It is an important physiological parameter involved in respiration and circulation. The Pulse Oximeter feature herein is small, portable, non-invasive and easy to use. The user only needs to insert a finger into the chamber to measure his/her SpO2 and pulse rate. The oximeter is suitable for home healthcare environments.
DEVICE INTENDED PURPOSE	The pulse oximeter is a reusable device and intended for spot-checking of oxygen saturation and pulse rate for use with the finger of adult patients.
VALIDITY DATE OF REGISTRATION	10/08/2026 - 09/08/2031

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
No results found.		