



DEVICE DETAILS

NAME OF DEVICE	INFRARED THERMOMETER
ESTABLISHMENT NAME	IHB HEALTHCARE SDN BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GB4696320-40017
BRAND NAME	SHENZHEN AEON TECHNOLOGY
MEDICAL DEVICE CATEGORY	MD 1100 - GENERAL ACTIVE MEDICAL DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	Shenzhen Aeon Technology Infrared Thermometer (A200 and A201) is an infrared thermometer, based on Infrared and digital technology, is intended for non-contact temperature detection. The product is suitable for use in family, hospital and clinic. Any object above the absolute zero (-273 ?) or more will be issued out of infrared, infrared thermometer through the sensor to receive the temperature data. Infrared thermometer used to measure human temperature reference design, using a very simple and convenient.
DEVICE INTENDED PURPOSE	Shenzhen Aeon Technology Infrared Thermometer (A200 and A201) is intended for the intermittent measurement and monitoring of human body temperature from forehead measurement at home, clinics and hospital.
VALIDITY DATE OF REGISTRATION	05/02/2020 - 04/02/2025

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Infrared Thermometer	A200
2	Infrared Thermometer	A201