



DEVICE DETAILS

NAME OF DEVICE	PULSE OXIMETER OLV-2700K
ESTABLISHMENT NAME	NIHON KOHDEN MALAYSIA SDN BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GC14352301217
BRAND NAME	NIHON KOHDEN (OLV-2700K)
MEDICAL DEVICE CATEGORY	MD 1302 - Monitoring devices of vital physiological parameters
DEVICE GROUPING TYPE	SYSTEM
DEVICE DESCRIPTION	THE OLV-2700K PULSE OXIMETER IS A PORTABLE INDOOR-USE MONITOR THAT MEASURES SPO2 OF A PATIENT IN THE HOSPITAL. THIS PULSE OXIMETER DISPLAYS THE PATIENT'S SPO2 AND PULSE RATE, AND GENERATES ALARMS. YOU CAN RUN IT ON AC OR BATTERY POWER ACCORDING TO THE LOCATION, SUCH AS ER, OPERATION ROOM, ICU, GENERAL WARD, ETC. THE PULSE OXIMETER IS DESIGNED SO THE OPERATOR CAN DIRECTLY PRESS THE KEYS FROM THE OPERATOR POSITION. THE LED DISPLAY SHOWS SPO2 VALUE AND PULSE RATE.
DEVICE INTENDED PURPOSE	THIS DEVICE MONITORS THE RESPIRATORY CONDITION OF A PATIENT BY MEASURING THE ARTERIAL OXYGEN SATURATION (SPO2) AND PULSE RATE (PR) IN A CONTINUOUS NON-INVASIVE MANNER. THE WAVEFORM SIGNALS CORRESPONDING TO THE VALUES MEASURED BY THIS DEVICE CAN BE ALSO OUTPUT BY CONNECTING TO A NIHON KOHDEN VITAL INFORMATION MONITOR VIA MONITOR INTERFACE KIT.
VALIDITY DATE OF REGISTRATION	25/08/2017 - 24/08/2022

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Oxypal Pulse Oximeter, OLV-2700K	OLV-2700K
2	Power cord H, 186656A	186656A
3	Power cord N, 314839B	314839B
4	Grounding lead, 544582A	544582A
5	Battery pack, YS-085P6	YS-085P6
6	Monitor interface kit, YS-085P7	YS-085P7

NO	NAME OF DEVICE	IDENTIFIER
7	Monitor interface kit, YS-085P8	YS-085P8
8	SpO2 connection cable, JL-302T	JL-302T
9	Finger probe BluPRO, TL-201T	TL-201T
10	Finger probe BluPRO, TL-630T1	TL-630T1
11	Finger probe BluPRO, TL-630T3	TL-630T3
12	Finger probe BluPRO, TL-631T1	TL-631T1
13	Finger probe BluPRO, TL-631T3	TL-631T3
14	Multi-site probe, TL-220T	TL-220T
15	SpO2 probe (for adult) BluPRO, TL-271T	TL-271T
16	SpO2 probe (for adult) BluPRO, TL-271T3	TL-271T3
17	SpO2 probe (for child) BluPRO, TL-272T	TL-272T
18	SpO2 probe (for child) BluPRO, TL-272T3	TL-272T3
19	SpO2 probe (for adult and neonate) BluPRO, TL-273T	TL-273T
20	SpO2 probe (for adult and neonate) BluPRO, TL-273T3	TL-273T3
21	SpO2 probe (for child and infant) BluPRO, TL-274T	TL-274T
22	SpO2 probe (for child and infant) BluPRO, TL-274T3	TL-274T3
23	Multi-site Y probe (disposable, for adult, child, neonate), TL-260T	TL-260T
24	Disposable SpO2 probe (for adult and neonate), TL-051S	TL-051S
25	Disposable SpO2 probe (for adult and neonate), TL-052S	TL-052S
26	Disposable SpO2 probe (for adult, child and infant), TL-061S	TL-061S
27	Disposable SpO2 probe (for adult, child and infant), TL-062S	TL-062S
28	Disposable SpO2 probe (for adult), TL-271T	TL-271T
29	Disposable SpO2 probe (for adult), TL-271T3	TL-271T3
30	Disposable SpO2 probe (for child), TL-272T	TL-272T
31	Disposable SpO2 probe (for child), TL-272T3	TL-272T3
32	Disposable SpO2 probe (for adult and neonate), TL-273T	TL-273T

NO	NAME OF DEVICE	IDENTIFIER
33	Disposable SpO2 probe (adult and neonate), TL-273T3	TL-273T3
34	Disposable SpO2 probe (for infant), TL-274T	TL-274T
35	Disposable SpO2 probe (for infant), TL-274T3	TL-274T3
36	Attachment tape, P263	P263
37	Foam tape, P260	P260
38	Cottony tape, 340703	P259
39	Sponge attachment tape S, P260A	P260A
40	Sponge attachment tape L, P260B	P260B
41	TL-260T clip adapter, P256	P256



MEDICAL DEVICE AUTHORITY