

DEVICE DETAILS			
NAME OF DEVICE	TEMP-PULSE OXIMETER		
ESTABLISHMENT NAME	AL GLOBAL CORPORATION (M) SDN BHD		
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE		
REGISTRATION NO	GC10725421-52611		
BRAND NAME	TEMP-PULSE OXIMETER/MEDLINKET		
MEDICAL DEVICE CATEGORY	MD 1100 - GENERAL ACTIVE MEDICAL DEVICES		
DEVICE GROUPING TYPE	SINGLE		
DEVICE DESCRIPTION	The Pulse Oximeter is intended for continuous use or spot checking in measuri and displaying functional arterial oxygen saturation (SpO2), pulse rate and temperature of patients in hospitals, physician's office, clinical settings and home care environment. The Temp-Pulse Oximeter can be used as a oxygen saturation monitor, heart rate monitor, body temperature thermometer and heal monitor for home daily use. And also can be a fitness trackers for sportspeople wishing to track their performance and fitness; Recognized Brand and Manufacturer by European Market: MED LINKET produce and export		
DEVICE INTENDED PURPOSE	Shenzhen Med-link Temp Pulse Oximeter is intended for continuous use or spot checking in measuring and displaying functional arterial oxygen saturation (SpO2),pulse rate and temperature of patients in hospitals,physician's office,clinical settings and home care environment. Target population:Adult,adolescent and child.		
VALIDITY DATE OF REGISTRATION	18/12/2020 - 17/12/2025		

LIST OF DEVICE			
NO	NAME OF DEVICE	IDENTIFIER	
1	TEMP-PULSE OXIMETER	AM-801	