



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	SOFT CONTACT LENS
<b>ESTABLISHMENT NAME</b>	MAXVUE VISION SDN BHD
<b>ROLE OF ESTABLISHMENT</b>	MANUFACTURER
<b>REGISTRATION NO</b>	GC172941046418
<b>BRAND NAME</b>	PHANTASEE
<b>MEDICAL DEVICE CATEGORY</b>	MD 0105 - Non-active ophthalmologic devices
<b>DEVICE GROUPING TYPE</b>	FAMILY
<b>DEVICE DESCRIPTION</b>	MAXVUE VISION SDN BHD'S SOFT CONTACT LENSES ARE HEMISPHERICAL SHELLS WITH MOULDED SPHERICAL BASE CURVES AND LATHE-CUT FRONT SURFACES. THE LENSES ARE FABRICATED FROM A NON-IONIC POLYMER. THE NON-IONIC LENS MATERIAL IS A HYDROPHILIC POLYMER OF 2-HYDROXYETHYL METHACRYLATE (2-HEMA) AND CROSS-LINKED WITH CROSSLINKING AGENT. THE LENSES ARE IMMERSUED IN NORMAL BUFFERED SALINE SOLUTION.
<b>DEVICE INTENDED PURPOSE</b>	MAXVUE VISIONS SOFT CONTACT LENSES FOR DAILY WEAR ARE INDICATED FOR THE CORRECTIONS OF VISUAL ACUITY IN APHAKIC AND NOT APHAKIC PERSON WITH NON-DISEASED EYES WITH MYOPIA OR HYPEROPIA. THE LENS MAY BE WORN BY PERSONS WHO EXHIBIT REFRACTIVE ASTIGMATISM OF 0.50 DIOPTERS OR LESS WHERE THE ASTIGMATISM DOES NOT INTERFERE WITH VISUAL ACUITY. THE LENS IS ALSO USED TO ENHANCE OR ALTER THE APPARENT COLOUR OF THE EYE. EYE CARE PRACTITIONERS MAY PRESCRIBE THE ABOVE LENSES FOR FREQUENT/PLANNED REPLACEMENT WEAR, WITH CLEANING DISINFECTION AND SCHEDULED REPLACEMENT. WHEN PRESCRIBED FOR FREQUENT/PLANNED REPLACEMENT WEAR, THE LENS MAY BE DISINFECTED USING A CHEMICAL DISINFECTING SYSTEM.
<b>VALIDITY DATE OF REGISTRATION</b>	02/07/2018 - 01/07/2023

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Phantasee Big Eyes Angel Hazel	9555644820787
2	Phantasee Big Eyes Brilliant Black	9555644820527
3	Phantasee Big Eyes Maya Blue	9555644823160
4	Phantasee Big Eyes Baby Brown	9555644820268

NO	NAME OF DEVICE	IDENTIFIER
5	Phantasee Big Eyes Pearl Grey	9555644820008
6	Phantasee Big Eyes Paris Green	9555644823429
7	Phantasee Big Eyes Beautiful Blue	9555644815677
8	Phantasee Big Eyes Charming Brown	9555644815417
9	Phantasee Big Eyes Defined Ring	9555644816452
10	Phantasee Big Eyes Passionate Purple	9555644815936
11	Phantasee Big Eyes Lovely Grey	9555644815158
12	Phantasee Big Eyes Lustrous Green	9555644814892
13	Phantasee Diva Blue	9555644818128
14	Phantasee Diva Choco	9555644818111
15	Phantasee Diva Grey	9555644818098
16	Phantasee Diva Green	9555644818081
17	Phantasee Diva Black	9555644818074
18	Phantasee Diva Brown	955564481810
19	Phantasee Natural Green	9555644818012
20	Phantasee Natural Hazel	9555644818043
21	Phantasee Natural Grey	9555644818029
22	Phantasee Natural Brown	9555644818036
23	Phantasee Natural Turquoise	9555644818067
24	Phantasee Natural Blue	9555644818050
25	Phantasee Vivid Blue	9555644817497
26	Phantasee Vivid Green	9555644816711
27	Phantasee Vivid Brown	9555644817237
28	Phantasee Vivid Grey	9555644816971
29	Phantasee Vivid Violet	9555644817756
30	Phantasee Fancy 14 MM	Bizarro

<b>NO</b>	<b>NAME OF DEVICE</b>	<b>IDENTIFIER</b>
31	Phantasee Fancy 14 MM	Blastaar
32	Phantasee Fancy 14 MM	Gaia
33	Phantasee Fancy 14 MM	Poseidon
34	Phantasee Fancy 14 MM	Neji Hyuga
35	Phantasee Fancy 14 MM	Rogue
36	Phantasee Fancy 14 MM	Snowbird
37	Phantasee Fancy 14 MM	Starfire
38	Phantasee Fancy 14 MM	Thunderbolts
39	Phantasee Fancy 14 MM	Autumn
40	Phantasee Fancy 14 MM	Black Cross
41	Phantasee Fancy 14 MM	Blue Artist
42	Phantasee Fancy 14 MM	Blue Mist
43	Phantasee Fancy 14 MM	Blue Sapphire
44	Phantasee Fancy 14 MM	Blue Tornado
45	Phantasee Fancy 14 MM	Blood Spin
46	Phantasee Fancy 14 MM	Charmed
47	Phantasee Fancy 14 MM	Chequered
48	Phantasee Fancy 14 MM	Dark Side
49	Phantasee Fancy 14 MM	Devil Cat
50	Phantasee Fancy 14 MM	Dragon Eye
51	Phantasee Fancy 14 MM	Dracula
52	Phantasee Fancy 14 MM	Green Galaxy
53	Phantasee Fancy 14 MM	Grim Reaper
54	Phantasee Fancy 14 MM	Kakashi
55	Phantasee Fancy 14 MM	Lestat Vampire
56	Phantasee Fancy 14 MM	Lizard Eye

<b>NO</b>	<b>NAME OF DEVICE</b>	<b>IDENTIFIER</b>
57	Phantasee Fancy 14 MM	Lost Galaxy
58	Phantasee Fancy 14 MM	Lunar Eclipse
59	Phantasee Fancy 14 MM	Green Mist
60	Phantasee Fancy 14 MM	Green Spin
61	Phantasee Fancy 14 MM	Green Tornado
62	Phantasee Fancy 14 MM	Green
63	Phantasee Fancy 14 MM	Madara
64	Phantasee Fancy 14 MM	Mangekyu
65	Phantasee Fancy 14 MM	Marijuana Leaf
66	Phantasee Fancy 14 MM	Pacific Blue
67	Phantasee Fancy 14 MM	Pentagram
68	Phantasee Fancy 14 MM	Pink Artist
69	Phantasee Fancy 14 MM	Pink Galaxy
70	Phantasee Fancy 14 MM	Pinky
71	Phantasee Fancy 14 MM	Pure Galaxy
72	Phantasee Fancy 14 MM	Purple
73	Phantasee Fancy 14 MM	Radiate
74	Phantasee Fancy 14 MM	Rage
75	Phantasee Fancy 14 MM	Red Galaxy
76	Phantasee Fancy 14 MM	Red Spin
77	Phantasee Fancy 14 MM	Red Tornado
78	Phantasee Fancy 14 MM	Rock Star
79	Phantasee Fancy 14 MM	Sasuke
80	Phantasee Fancy 14 MM	Saw
81	Phantasee Fancy 14 MM	Snake Eye
82	Phantasee Fancy 14 MM	Soccer

<b>NO</b>	<b>NAME OF DEVICE</b>	<b>IDENTIFIER</b>
83	Phantasee Fancy 14 MM	Spider
84	Phantasee Fancy 14 MM	Summer
85	Phantasee Fancy 14 MM	Temptress
86	Phantasee Fancy 14 MM	Twilight
87	Phantasee Fancy 14 MM	Volcano
88	Phantasee Fancy 14 MM	Volturi
89	Phantasee Fancy 14 MM	Werewolf
90	Phantasee Fancy 14 MM	White Galaxy
91	Phantasee Fancy 14 MM	White Manson
92	Phantasee Fancy 14 MM	White Tornado
93	Phantasee Fancy 14 MM	White Witch
94	Phantasee Fancy 14 MM	Winter
95	Phantasee Fancy 14 MM	Yellow Artist
96	Phantasee Fancy 14 MM	Yellow Cat
97	Phantasee Fancy 14 MM	Yellow Galaxy
98	Phantasee Fancy 14 MM	Yellow Spin
99	Phantasee Fancy 14 MM	Yellow Tornado
100	Phantasee Fancy 14 MM	Yellow
101	Phantasee Fancy 17 MM	Black Titan
102	Phantasee Fancy 17 MM	Daredevil
103	Phantasee Fancy 22 MM	Sabretooth
104	Phantasee Fancy 22 MM	Cyclop
105	Phantasee Fancy 22 MM	Medusa
106	Phantasee Fancy 22 MM	Gremlin
107	Phantasee Fancy 22 MM	Athena Blue
108	Phantasee Fancy 22 MM	Caliban

NO	NAME OF DEVICE	IDENTIFIER
109	Phantasee Fancy 22 MM	Raiden Pink



**MEDICAL DEVICE AUTHORITY**