




DEVICE DETAILS

NAME OF DEVICE	DISPOSABLE PROTECTIVE SURGICAL FACE MASK
ESTABLISHMENT NAME	EMPEROR DYNASTY CORPORATION SDN BHD
ROLE OF ESTABLISHMENT	MANUFACTURER
REGISTRATION NO	GA8138620-50273
BRAND NAME	MEDICURVE
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	3-PLY EARLOOP SURGICAL FACE MASK; LATEX FREE; FLUID RESISTANT; BACTERIAL FILTRATION EFFICIENCY > 99%, PARTIAL FILTRATION EFFICIENCY >95%
DEVICE INTENDED PURPOSE	SURGICAL FACE MASK TO PREVENT DISEASE THROUGH FLUID AND AIR TRANSMITTANCE, NON-INVASIVE AND NON-ACTIVE MEDICAL DEVICE
VALIDITY DATE OF REGISTRATION	03/12/2020 - 02/12/2025

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	DISPOSABLE PROTECTIVE SURGICAL FACE MASK 	MD10001
2	DISPOSABLE PROTECTIVE SURGICAL FACE MASK	MD10002