




DEVICE DETAILS

NAME OF DEVICE	DISPOSABLE SURGICAL FACE MASK
ESTABLISHMENT NAME	SIMPLY K SDN BHD
ROLE OF ESTABLISHMENT	MANUFACTURER
REGISTRATION NO	GA4306820-50594
BRAND NAME	MEDICMAD
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	Disposable surgical face masks embedded with nose clip and ear loops or with ties. This surgical face mask made of 3-ply of premium non-woven fabric with Bacterial Filtration Efficiency (BFE) >99% and Particle Filtration Efficiency (PFE) >95%.
DEVICE INTENDED PURPOSE	Disposable Surgical face mask that covers the user's nose and mouth and provides a physical barrier from transmission of viruses, microorganisms and body fluids. It is suitable to be use for surgical procedures to protect both patient and healthcare provider.
VALIDITY DATE OF REGISTRATION	16/12/2020 - 15/12/2025

LIST OF DEVICE

NO	NAME OF DEVICE 	IDENTIFIER
1	DISPOSABLE ADULT SURGICAL FACE MASK	ADULT SURGICAL FACE MASK
2	DISPOSABLE CHILDREN SURGICAL FACE MASK	CHILDREN SURGICAL FACE MASK