



DEVICE DETAILS

NAME OF DEVICE	3 PLY SURGICAL FACE MASK
ESTABLISHMENT NAME	PINANG MEDICAL SUPPLIES SDN. BHD.
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GA4807920-51157
BRAND NAME	LIFEPLUS
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	A surgical mask is a loose-fitting, disposable device that creates a physical barrier between the mouth and nose of the wearer and potential contaminants in the immediate environment
DEVICE INTENDED PURPOSE	Fluid resistant and provides the wearer protection against large droplets, splashes, or sprays of bodily or other hazardous fluids. Protects the patient from the wearer's respiratory emissions.
VALIDITY DATE OF REGISTRATION	29/12/2020 - 28/12/2025

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	SURGICAL MASK 3-PLY (50'S): LIFEPLUS (EAR-LOOP)	EAR-LOOP
2	SURGICAL MASK 3-PLY (50'S): LIFEPLUS (TIE-ON)	TIE-ON
3	SURGICAL MASK 3-PLY (50'S): LIFEPLUS (HIJAB)	HIJAB
4	SURGICAL MASK 3-PLY (50'S): LIFEPLUS (KIDS) (EAR-LOOP)	EAR-LOOP