

DEVICE DETAILS		
NAME OF DEVICE	FINGERTIP PULSE OXIMETER	
<b>ESTABLISHMENT NAME</b>	GAS WALKER MIRAMED SDN BHD	
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE	
REGISTRATION NO	GC9480821-57718	
BRAND NAME	YUWELL	
MEDICAL DEVICE CATEGORY	MD 1100 - GENERAL ACTIVE MEDICAL DEVICES	
<b>DEVICE GROUPING TYPE</b>	FAMILY	
DEVICE DESCRIPTION	Pulse oximeter is a noninvasive method for monitoring a person's oxygen saturation (SO2). Though its reading of SpO2 (peripheral oxygen saturation) is not always identical to the more desirable reading of SaO2 (arterial oxygen saturation) from arterial blood gas analysis, the two are correlated well enough that the safe, convenient, noninvasive, inexpensive pulse oximetry method is valuable for measuring oxygen saturation in clinical use.	
DEVICE INTENDED PURPOSE	The fingertip Oximeter can be used to measure human Haemoglobin Saturation and heart rate through finger, the product is suitable for family, hospital (including clinical use in internist/surgery, Anaesthesia, paediatrics, intensive sports (It can be used before or after sport. Operation in sport is not recommended) and etc. The product is not suitable to monitor patient continuously.	
VALIDITY DATE OF REGISTRATION	24/03/2021 - 23/03/2026	

LIST OF DEVICE			
NO	NAME OF DEVICE	IDENTIFIER	
1	Fingertip Pulse Oximeter	YX302	
2	Fingertip Pulse Oximeter	YX100	