



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	XP SERIES
<b>ESTABLISHMENT NAME</b>	VENUSYS MEDICAL SDN BHD
<b>ROLE OF ESTABLISHMENT</b>	AUTHORIZED REPRESENTATIVE
<b>REGISTRATION NO</b>	GC8237523616
<b>BRAND NAME</b>	FOTONA
<b>MEDICAL DEVICE CATEGORY</b>	MD 1202 - Imaging devices utilising non-ionizing radiation
<b>DEVICE GROUPING TYPE</b>	FAMILY
<b>DEVICE DESCRIPTION</b>	THIS DEVICE IS A MEDICAL SURGICAL LASER OR AESTHETIC LASER CONSIST OF A WAVELENGTH 1064NM, A LONG-PULSED ND:YAG LASER WHICH IS MEANT FOR APPLICATIONS IN DERMATOLOGY AND GYNECOLOGY.
<b>DEVICE INTENDED PURPOSE</b>	THIS MEDICAL DEVICE IS USED IN DERMATOLOGY AND GYNECOLOGY TREATMENT SUCH AS VASCULAR LESIONS AND VAGINAL TIGHTENING.
<b>VALIDITY DATE OF REGISTRATION</b>	14/05/2021 - 13/05/2026

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	FOTONA XP DYNAMIS	85992 & 92730 (Touch Screen)
2	FOTONA XP SPECTRO	87005 & 92732 (Touch Screen)