



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	Disposable Medical Face Mask
<b>ESTABLISHMENT NAME</b>	FERMETEC MEDICAL SDN BHD
<b>ROLE OF ESTABLISHMENT</b>	AUTHORIZED REPRESENTATIVE
<b>REGISTRATION NO</b>	GA4445721-62014
<b>BRAND NAME</b>	Hubei Kangning
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	FAMILY
<b>DEVICE DESCRIPTION</b>	The product consists of a mask body clip and a mask belt. The mask body is three-layer structure, the inner and outer layers are made of polypropylene (PP) non-woven fabric, the intermediate layer filter material is polypropylene fused spray cloth clip made of metal coated plastic material, and the mask belt is made of elastic material.
<b>DEVICE INTENDED PURPOSE</b>	The Medical face mask are intended to be worn it protect against the spread of transmission of infections germ patient against infectious germs. In addition, in certain situations the wearer should be protected against splashes of the potentially contaminated liquids and viable particles.
<b>VALIDITY DATE OF REGISTRATION</b>	24/05/2021 - 23/05/2026

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Medical Mask	KN001-01
2	Medical Mask	KN001-02