



## DEVICE DETAILS

NAME OF DEVICE	Disposable Medical Face Mask
ESTABLISHMENT NAME	FERMETEC MEDICAL SDN BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GA4445721-62014
BRAND NAME	Hubei Kangning
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	The product consists of a mask body clip and a mask belt. The mask body is three-layer structure, the inner and outer layers are made of polypropylene (PP) non-woven fabric, the intermediate layer filter material is polypropylene fused spray cloth clip made of metal coated plastic material, and the mask belt is made of elastic material.
DEVICE INTENDED PURPOSE	The Medical face mask are intended to be worn it protect against the spread of transmission of infections germ patient against infectious germs. In addition, in certain situations the wearer should be protected against splashes of the potentially contaminated liquids and viable particles.
VALIDITY DATE OF REGISTRATION	24/05/2021 - 23/05/2026

## LIST OF DEVICE



NO	NAME OF DEVICE	IDENTIFIER
1	Medical Mask	KN001-01
2	Medical Mask	KN001-02