



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	3PLY Disposable Surgical Face mask
<b>ESTABLISHMENT NAME</b>	PEN TELL MARKETING SDN BHD
<b>ROLE OF ESTABLISHMENT</b>	AUTHORIZED REPRESENTATIVE
<b>REGISTRATION NO</b>	GA9366521-63698
<b>BRAND NAME</b>	ZHUSHI
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	SINGLE
<b>DEVICE DESCRIPTION</b>	this product composed of 3 materials: non-woven fabric, nose strip and elastic band. the face mask devided into inner, middle and outer layers,Blue Colour , Size 17.5cm x 9.5cm
<b>DEVICE INTENDED PURPOSE</b>	For filtering particulate matter in the air, blocking droplets, blood, body fluids and secretions in the medical working environment, and shall not be used to isolate rooms in intensive care units and other places with strict requirements for the control of microbial indicators.
<b>VALIDITY DATE OF REGISTRATION</b>	15/06/2021 - 14/06/2026

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	3PLY Disposable Surgical Face mask	ZS-B