



DEVICE DETAILS

| | |
|--------------------------------------|---|
| NAME OF DEVICE | DISPOSABLE MEDICAL FACE MASK |
| ESTABLISHMENT NAME | SIMPLY K SDN BHD |
| ROLE OF ESTABLISHMENT | MANUFACTURER |
| REGISTRATION NO | GA10611221-65198 |
| BRAND NAME | EHOP |
| MEDICAL DEVICE CATEGORY | SINGLE-USE DEVICES |
| DEVICE GROUPING TYPE | FAMILY |
| DEVICE DESCRIPTION | Disposable medical face masks embedded with nose clip and ear loops or with head loops. This medical face mask made of 3-ply of premium non-woven fabric with Bacterial Filtration Efficiency (BFE) >95% and Particle Filtration Efficiency (PFE) >95%. |
| DEVICE INTENDED PURPOSE | Disposable Medical face mask that covers the user's nose and mouth and provides a physical barrier from transmission of viruses, microorganisms and dust. It is intended to be use by healthcare provider and other persons for medical purposes and to control infectious diseases, particularly in epidemic or pandemic situations. |
| VALIDITY DATE OF REGISTRATION | 28/06/2021 - 27/06/2026 |

LIST OF DEVICE

| NO | NAME OF DEVICE | IDENTIFIER |
|----|---------------------------------------|--|
| 1 | DISPOSABLE ADULT MEDICAL FACE MASK | ADULT MEDICAL FACE MASK - EAR LOOP |
| 2 | DISPOSABLE ADULT MEDICAL FACE MASK | ADULT MEDICAL FACE MASK - HEAD LOOP |
| 3 | DISPOSABLE CHILDREN MEDICAL FACE MASK | CHILDREN MEDICAL FACE MASK - EAR LOOP |
| 4 | DISPOSABLE CHILDREN MEDICAL FACE MASK | CHILDERN MEDICAL FACE MASK - HEAD LOOP |