




## DEVICE DETAILS

NAME OF DEVICE	Non Woven Surgical Face Mask (Ear Loop)
ESTABLISHMENT NAME	QUANTUM UPSTREAM SDN.BHD.
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GA3775521-65292
BRAND NAME	HospiGuard
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	Non Woven Surgical Face Mask (Ear Loop) is made from 2 layers of non-woven fabric and 1 layer of melt-blown filter fabric with high BFE over 95%, hypoallergenic and very low resistance to breath
DEVICE INTENDED PURPOSE	Surgical Face Mask is intended for used as a barrier to prevent air-borne transfer of bacteria and particles. This is intended for used in operating theater.
VALIDITY DATE OF REGISTRATION	28/06/2021 - 27/06/2026

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Non Woven Surgical Face Mask (Ear Loop) 	BP01
2	Non Woven Surgical Face Mask (Ear Loop)	MP01
3	Non Woven Surgical Face Mask (Ear Loop)	FL01
4	Non Woven Surgical Face Mask (Ear Loop)	MG01
5	Non Woven Surgical Face Mask (Ear Loop)	OA10
6	Non Woven Surgical Face Mask (Ear Loop)	CL01
7	Non Woven Surgical Face Mask (Ear Loop)	CV01
8	Non Woven Surgical Face Mask (Ear Loop)	NB01