



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	SURGICAL FACE MASK DISPOSABLE
<b>ESTABLISHMENT NAME</b>	EUGIN RESOURCES SDN BHD
<b>ROLE OF ESTABLISHMENT</b>	MANUFACTURER
<b>REGISTRATION NO</b>	GA3934521-65298
<b>BRAND NAME</b>	HENNSON
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	SINGLE
<b>DEVICE DESCRIPTION</b>	This product consists of a mask body (bottom layer, top layer and middle layer), a nose clip and a mask band. The bottom and top layer of the mask body are made of breathable, soft and skin friendly non-woven fabric. The middle layer is made of polypropylene melt blown cloths. The mask is equipped with a nose clip which made of a pliable metal strip. The mask belt is made of polyester and non-woven fabric. It is suitable for all types of clinical personnel to wear during non-invasive procedures, covering the user's mouth, nose and jaw.
<b>DEVICE INTENDED PURPOSE</b>	Surgical mask is intended for single use by operating room personnel and other general healthcare workers to protect both patients and healthcare workers against transfer of microorganisms, blood and body fluids, and airborne particulates.
<b>VALIDITY DATE OF REGISTRATION</b>	28/06/2021 - 27/06/2026

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	SURGICAL FACE MASK DISPOSABLE	SFMD