



DEVICE DETAILS

NAME OF DEVICE	PULSE OXIMETER
ESTABLISHMENT NAME	PERNIAGAAN BEKALAN ABADI SDN BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GB9253021-66450
BRAND NAME	LITEXIM
MEDICAL DEVICE CATEGORY	MD 1100 - GENERAL ACTIVE MEDICAL DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	The Fingertip Pulse Oximeter is a non-invasive device intended for the spot-check of oxygen saturation of arterial hemoglobin (SpO ₂) and the pulse rate of adult and testes in home This device is not intended for continuous monitoring. The pulse oxygen saturation is the percentage of HbO ₂ in the total Hb in the blood, so called the O ₂ concentration in the blood. It is an important bic-parameter for the respiration. For the purpose of measuring the SpO ₂ more and accurately, our company developed the Pulse Oximeter. At the same time, the device can measure the pulse rate simultaneously. The Pulse Oximeter features in small volume, low power consumption, convenient operation and being portable. It is only necessary for testes to put one of his fingers into a fingertip photoelectric sensor for diagnosis, and a display screen will directly show measured value of Hemoglobin Saturation.
DEVICE INTENDED PURPOSE	The Fingertip Pulse Oximeter is a non-invasive device intended for the spot-check of oxygen saturation of arterial hemoglobin (SpO ₂) and the pulse rate of adult and testes in home This device is not intended for continuous monitoring
VALIDITY DATE OF REGISTRATION	21/06/2021 - 20/06/2026

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Pulse Oximeter	OX206
2	Pulse Oximeter	OX205