

DEVICE DETAILS		
NAME OF DEVICE	PULSE OXIMETER	
ESTABLISHMENT NAME	PERNIAGAAN BEKALAN ABADI SDN BHD	
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE	
REGISTRATION NO	GB9253021-66450	
BRAND NAME	LITEXIM	
MEDICAL DEVICE CATEGORY	MD 1100 - GENERAL ACTIVE MEDICAL DEVICES	
DEVICE GROUPING TYPE	FAMILY	
DEVICE DESCRIPTION	The Fingertip Pulse Oximeter is a non-invasive device intended for the spotcheck of oxygen saturation of arterial hemoglobin (SpO2) and the pulse rate of adult and testes in home This device is not intended for continuous monitoring. The pulse oxygen saturation is the percentage of HbO2 in the total Hb in the blood, so called the OZ concentration in the blood. It is an important bic-parameter for the respiration. For the purpose of measuring the SpO2 more and accurately, our company developed the Pulse Oximeter. At the same time, the device can measure the pulse rate simultaneously. The Pulse Oximeter features in small volume, low power consumption, convenient operation and being portable. It is only necessary for testes to put one of his fingers into a fingertip photoelectric sensor for diagnosis, and a display screen will directly show measured value of Hemoglobin Saturation.	
DEVICE INTENDED PURPOSE	The Fingertip Pulse Oximeter is a non-invasive device intended for the spot- check of oxygen saturation of arterial hemoglobin (SpO2) and the pulse rate of adult and testes in home This device is not intended for continuous monitoring	
VALIDITY DATE OF REGISTRATION	21/06/2021 - 20/06/2026	

LIST OF DEVICE			
NO	NAME OF DEVICE	IDENTIFIER	
1	Pulse Oximeter	OX206	
2	Pulse Oximeter	OX205	