



DEVICE DETAILS

NAME OF DEVICE	FINGER PULSE OXIMETER
ESTABLISHMENT NAME	PFLEGE MEDICAL SUPPLIES SDN. BHD.
ROLE OF ESTABLISHMENT	MANUFACTURER
REGISTRATION NO	GC5747021-66902
BRAND NAME	PROFIX
MEDICAL DEVICE CATEGORY	MD 1302 - Monitoring devices of vital physiological parameters
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	PROFIX Finger Pulse Oximeter are fingertip pulse oximeter, can be used to measure human Hemoglobin Saturation and pulse rate through finger. The product is suitable for family, hospital (including clinical use in internist/surgery, anesthesia, pediatrics, intensive care and etc.), oxygen bar, social medical organizations and physical care in sports.
DEVICE INTENDED PURPOSE	PROFIX Finger Pulse Oximeter is intended for non-invasive monitoring of the adult or child's SpO2 and PR.
VALIDITY DATE OF REGISTRATION	17/07/2021 - 16/07/2026

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Finger Pulse Oximeter	M130
2	Finger Pulse Oximeter	M160
3	Finger Pulse Oximeter	M170