



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	4-PLY SURGICAL NON-STERILE DISPOSABLE FACE MASK
<b>ESTABLISHMENT NAME</b>	EUGIN RESOURCES SDN BHD
<b>ROLE OF ESTABLISHMENT</b>	MANUFACTURER
<b>REGISTRATION NO</b>	GA8165521-66888
<b>BRAND NAME</b>	EUGINMask
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	SINGLE
<b>DEVICE DESCRIPTION</b>	Surgical mask is a single use 4 layers mask with outer and inner layer (spunbond polypropylene) that sandwich two layers of melt-blown polypropylene filter material. Earloops are of urethane elastic fiber and latex-free. The nose bridge is a pliable white aluminium stripe, covered by PP covering. All of the materials used in the construction of the Surgical Mask are being used in current marketed device.
<b>DEVICE INTENDED PURPOSE</b>	Surgical mask is intended for single use by operating room personnel and other general healthcare workers to protect both patients and healthcare workers against transfer of microorganisms, blood and body fluids, and airborne particulates.
<b>VALIDITY DATE OF REGISTRATION</b>	17/07/2021 - 16/07/2026

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	4-PLY SURGICAL NON-STERILE DISPOSABLE FACE MASK	4PSNSDFM