



DEVICE DETAILS

NAME OF DEVICE	4-PLY SURGICAL NON-STERILE DISPOSABLE FACE MASK
ESTABLISHMENT NAME	EUGIN RESOURCES SDN BHD
ROLE OF ESTABLISHMENT	MANUFACTURER
REGISTRATION NO	GA8165521-66888
BRAND NAME	EUGINMask
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	SINGLE
DEVICE DESCRIPTION	Surgical mask is a single use 4 layers mask with outer and inner layer (spunbond polypropylene) that sandwich two layers of melt-blown polypropylene filter material. Earloops are of urethane elastic fiber and latex-free. The nose bridge is a pliable white aluminium stripe, covered by PP covering. All of the materials used in the construction of the Surgical Mask are being used in current marketed device.
DEVICE INTENDED PURPOSE	Surgical mask is intended for single use by operating room personnel and other general healthcare workers to protect both patients and healthcare workers against transfer of microorganisms, blood and body fluids, and airborne particulates.
VALIDITY DATE OF REGISTRATION	17/07/2021 - 16/07/2026

LIST OF DEVICE



NO	NAME OF DEVICE	IDENTIFIER
1	4-PLY SURGICAL NON-STERILE DISPOSABLE FACE MASK	4PSNSDFM