



## DEVICE DETAILS

NAME OF DEVICE	HYGIENIC FACE MASKS (ADULTS & KIDS)
ESTABLISHMENT NAME	WATSON'S PERSONAL CARE STORES SDN BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GA10740721-67940
BRAND NAME	WATSONS
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	Face mask made of 3-ply pleated, non-woven material.
DEVICE INTENDED PURPOSE	Able to filter droplets, pollen, dusts and other air particulates. The filter layer of the face mask has a Bacterial Filtration Efficiency (BFE), Particle Filtration Efficiency (PFE) and Viral Filtration Efficiency (VFE) of more than 99% thus, providing the wearer protection from the transfer of particles from the environment.
VALIDITY DATE OF REGISTRATION	29/07/2021 - 28/07/2026

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Watsons Hygienic Face Mask (Adults)	65030
2	Watsons Hygienic Face Mask (Kids)	65031