



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	PULSE OXIMETER
<b>ESTABLISHMENT NAME</b>	QUANTUM UPSTREAM SDN.BHD.
<b>ROLE OF ESTABLISHMENT</b>	AUTHORIZED REPRESENTATIVE
<b>REGISTRATION NO</b>	GC10889321-68404
<b>BRAND NAME</b>	EZYSWISS
<b>MEDICAL DEVICE CATEGORY</b>	MD 1302 - Monitoring devices of vital physiological parameters
<b>DEVICE GROUPING TYPE</b>	SINGLE
<b>DEVICE DESCRIPTION</b>	ezyswiss Pulse Oximeter is intended for non-invasive spot-check measurement of functional oxygen saturation of arterial hemoglobin (SpO <sub>2</sub> ). Advanced DSP algorithm can minimize the influence of motion artifact and improve measurement accuracy of low perfusion. The Pulse Oximeter can be used to measure human SpO <sub>2</sub> and heart rate through finger. The product is suitable for family, hospital (including clinical use in internist/surgery, Anesthesia, pediatrics etc.) Oxygen Bar, social medical organizations, physical care in sports etc.
<b>DEVICE INTENDED PURPOSE</b>	ezyswiss Pulse Oximeter Pulse Oximeter is a non-invasive device intended for spot checking of functional oxygen saturation of arterial hemoglobin (SpO <sub>2</sub> ) and pulse rate (PR). This portable device is indicated for use in adult patients in clinical institution and home environments.
<b>VALIDITY DATE OF REGISTRATION</b>	03/08/2021 - 02/08/2026

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	PULSE OXIMETER	OXI-01