



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	PULSE OXIMETER
<b>ESTABLISHMENT NAME</b>	PINANG MEDICAL SUPPLIES SDN. BHD.
<b>ROLE OF ESTABLISHMENT</b>	AUTHORIZED REPRESENTATIVE
<b>REGISTRATION NO</b>	GC5944221-70000
<b>BRAND NAME</b>	LIFEPLUS
<b>MEDICAL DEVICE CATEGORY</b>	MD 1302 - Monitoring devices of vital physiological parameters
<b>DEVICE GROUPING TYPE</b>	FAMILY
<b>DEVICE DESCRIPTION</b>	Lifepus Pulse Oximeter is a portable, convenient, non-invasive device, used for monitoring arterial haemoglobin oxygen saturation (SP02) and pulse rate. The personal application are adult patients (weight: >30kgs) and pediatric patients(weight:20-30kgs). Index finger, middle finger and ring finger are suitable position for monitor. It's intended for spot-checking or attended-care monitoring in Home Health Care and Medical Facility.
<b>DEVICE INTENDED PURPOSE</b>	Lifepus Pulse Oximeter is a portable, convenient, non-invasive device, intended for used for monitoring arterial haemoglobin oxygen saturation (SP02) and pulse rate. It's intended for spot-checking or attended-care monitoring in Home Health Care and Medical Facility.
<b>VALIDITY DATE OF REGISTRATION</b>	25/08/2021 - 24/08/2026

## LIST OF DEVICE

MEDICAL DEVICE AUTHORITY

NO	NAME OF DEVICE	IDENTIFIER
1	Pulse Oximeter	FS10K
2	Pulse Oximeter	FS20K
3	Pulse Oximeter	FS10I