



DEVICE DETAILS

NAME OF DEVICE	EXACTAMED Oral Dispenser
ESTABLISHMENT NAME	BAXTER HEALTHCARE (MALAYSIA) SDN. BHD.
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GMD8611690417A
BRAND NAME	EXACTAMED
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	Three-piece plastic dispenser of modified syringe design. Available in amber or clear. Each dispenser clearly marked with graduations.
DEVICE INTENDED PURPOSE	Exacta-Med Oral Liquid Dispensers are intended to be used to provide an accurate method of measuring and dispensing liquid medications to a patient orally. The Pharmacy Pack includes Oral Dispensers and Self-Righting Tip Caps. The dispenser is a three part syringe-like device designed to provide an accurate method of measuring and dispensing oral liquid medications. The dispenser tip is designed to not accept a standard luer taper medical needle. The dispenser provides superior volume accuracy for oral liquid medication doses versus traditional household teaspoon or other commonly used dosing devices
VALIDITY DATE OF REGISTRATION	16/01/2022 - 15/01/2027

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	EXACTANED Oral Dispensers, 10ml, Clear, Pharmacy Pack	H9387110
2	EXACTANED Oral Dispensers, 20ml, Clear, Pharmacy Pack	H9387120
3	EXACTANED Oral Dispensers, 60ml, Clear, Pharmacy Pack	H9387160