



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	EXACTAMED Oral Dispenser
<b>ESTABLISHMENT NAME</b>	BAXTER HEALTHCARE (MALAYSIA) SDN. BHD.
<b>ROLE OF ESTABLISHMENT</b>	AUTHORIZED REPRESENTATIVE
<b>REGISTRATION NO</b>	GMD8611690417A
<b>BRAND NAME</b>	EXACTAMED
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	FAMILY
<b>DEVICE DESCRIPTION</b>	Three-piece plastic dispenser of modified syringe design. Available in amber or clear. Each dispenser clearly marked with graduations.
<b>DEVICE INTENDED PURPOSE</b>	Exacta-Med Oral Liquid Dispensers are intended to be used to provide an accurate method of measuring and dispensing liquid medications to a patient orally. The Pharmacy Pack includes Oral Dispensers and Self-Righting Tip Caps. The dispenser is a three part syringe-like device designed to provide an accurate method of measuring and dispensing oral liquid medications. The dispenser tip is designed to not accept a standard luer taper medical needle. The dispenser provides superior volume accuracy for oral liquid medication doses versus traditional household teaspoon or other commonly used dosing devices
<b>VALIDITY DATE OF REGISTRATION</b>	16/01/2022 - 15/01/2027

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	EXACTANED Oral Dispensers, 10ml, Clear, Pharmacy Pack	H9387110
2	EXACTANED Oral Dispensers, 20ml, Clear, Pharmacy Pack	H9387120
3	EXACTANED Oral Dispensers, 60ml, Clear, Pharmacy Pack	H9387160