



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	3way stopcock with/without extension tube
<b>ESTABLISHMENT NAME</b>	COMCORDE MEDICAL (M) SDN BHD
<b>ROLE OF ESTABLISHMENT</b>	MANUFACTURER
<b>REGISTRATION NO</b>	GA3339619-26780
<b>BRAND NAME</b>	Hilfsmittel
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	FAMILY
<b>DEVICE DESCRIPTION</b>	Stopcocks are directional valves used to regulate flow of fluid through attached tubes. The stopcocks connections are 6% taper luer conical connectors at all ends of the device. Two of which are female Luer taper connections while the other one is a Male Luer taper connection with a rotating locking ring. The Stopcock comes with an extended variant where the Male connection has an additional flexible hollow tube attached to allow better flexibility of placing the medical device.
<b>DEVICE INTENDED PURPOSE</b>	The stopcocks is intended to be attached onto infusion lines for the purpose of regulating flow in an infusion line. Stopcock with extension line provides additional flexibility of placing the stopcock.
<b>VALIDITY DATE OF REGISTRATION</b>	11/03/2019 - 10/03/2024

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	3 way Stopcock (Blue)	CM-B3
2	3 way Stopcock (Red)	CM-R3
3	3 Way Stopcock (White)	CM-W3
4	3 Way Stopcock (Yellow)	CM-Y3
5	2mm Extension Tube 5cm with 3 Way stopcock Blue	CM-B3-ET02-05
6	2mm Extension Tube 10cm with 3 way stopcock blue	CM-B3-ET02-10
7	2mm Extension Tube 15cm with 3 way stopcock blue	CM-B3-ET02-15
8	2mm Extension Tube 30cm with 3 way stopcock blue	CM-B3-ET02-30
9	2mm Extension Tube 140cm with 3 way stopcock blue	CM-B3-ET02-140
10	2mm Extension Tube (blue line) 5cm with 3 way stopcock white	CM-B3-ET02-05B

NO	NAME OF DEVICE	IDENTIFIER
11	2mm Extension Tube (blue line) 10cm with 3 way stopcock blue	CM-B3-ET02-10B
12	2mm Extension Tube (blue line) 15cm with 3 way stopcock blue	CM-B3-ET02-15B
13	2mm Extension Tube (blue line) 30cm with 3 way stopcock blue	CM-B3-ET02-30B
14	2mm Extension Tube (blue line) 140cm with 3 way stopcock blue	CM-B3-ET02-140B
15	2mm Extension Tube 5cm with 3 way stopcock white	CM-W3-ET02-05
16	2mm Extension Tube 10cm with 3 way stopcock white	CM-W3-ET02-10
17	2mm Extension Tube 15cm with 3 way stopcock white	CM-W3-ET02-15
18	2mm Extension Tube 30cm with 3 way stopcock white	CM-W3-ET02-30
19	2mm Extension Tube 140cm with 3 way stopcock white	CM-W3-ET02-140
20	3mm Extension Tube 5cm with 3 Way stopcock Blue	CM-B3-ET03-05
21	3mm Extension Tube 10cm with 3 way stopcock blue	CM-B3-ET03-10
22	3mm Extension Tube 15cm with 3 way stopcock blue	CM-B3-ET03-15
23	3mm Extension Tube 30cm with 3 way stopcock blue	CM-B3-ET03-30
24	3mm Extension Tube 140cm with 3 way stopcock blue	CM-B3-ET03-140
25	3mm Extension Tube (blue line) 5cm with 3 Way stopcock Blue	CM-B3-ET03-05B
26	3mm Extension Tube (blue line) 10cm with 3 way stopcock blue	CM-B3-ET03-10B
27	3mm Extension Tube (blue line) 15cm with 3 way stopcock blue	CM-B3-ET03-15B
28	3mm Extension Tube (blue line) 30cm with 3 way stopcock blue	CM-B3-ET03-30B
29	3mm Extension Tube (blue line) 140cm with 3 way stopcock blue	CM-B3-ET03-140B
30	3mm Extension Tube 5cm with 3 way stopcock white	CM-W3-ET03-05
31	3mm Extension Tube 10cm with 3 way stopcock white	CM-W3-ET03-10
32	3mm Extension Tube 15cm with 3 way stopcock white	CM-W3-ET03-15
33	3mm Extension Tube 30cm with 3 way stopcock white	CM-W3-ET03-30
34	3mm Extension Tube 140cm with 3 way stopcock white	CM-W3-ET03-140
35	4mm Extension Tube 5cm with 3 Way stopcock Blue	CM-B3-ET04-05
36	4mm Extension Tube 10cm with 3 way stopcock blue	CM-B3-ET04-10

NO	NAME OF DEVICE	IDENTIFIER
37	4mm Extension Tube 15cm with 3 way stopcock blue	CM-B3-ET04-15
38	4mm Extension Tube 30cm with 3 way stopcock blue	CM-B3-ET04-30
39	4mm Extension Tube 140cm with 3 way stopcock blue	CM-B3-ET04-140
40	4mm Extension Tube (blue line) 5cm with 3 Way stopcock Blue	CM-B3-ET04-05B
41	4mm Extension Tube (blue line) 10cm with 3 way stopcock blue	CM-B3-ET04-10B
42	4mm Extension Tube (blue line) 15cm with 3 way stopcock blue	CM-B3-ET04-15B
43	4mm Extension Tube (blue line) 30cm with 3 way stopcock blue	CM-B3-ET04-30B
44	4mm Extension Tube (blue line) 140cm with 3 way stopcock blue	CM-B3-ET04-140B
45	4mm Extension Tube 5cm with 3 way stopcock white	CM-W3-ET04-05
46	4mm Extension Tube 10cm with 3 way stopcock white	CM-W3-ET04-10
47	4mm Extension Tube 15cm with 3 way stopcock white	CM-W3-ET04-15
48	4mm Extension Tube 30cm with 3 way stopcock white	CM-W3-ET04-30
49	4mm Extension Tube 140cm with 3 way stopcock white	CM-W3-ET04-140
50	4mm Extension Tube (white line) 5cm with 3 way stopcock white	CM-W3-ET04-05W
51	4mm Extension Tube (white line) 10cm with 3 way stopcock white	CM-W3-ET04-10W
52	4mm Extension Tube (white line) 15cm with 3 way stopcock white	CM-W3-ET04-15W
53	4mm Extension Tube (white line) 30cm with 3 way stopcock white	CM-W3-ET04-30W
54	4mm Extension Tube (white line) 140cm with 3 way stopcock white	CM-W3-ET04-140W
55	2mm extension pressure tube 5cm with 3 way stopcock red	CM-R3-PT02-05
56	2mm extension pressure tube 10cm with 3 way stopcock red	CM-R3-PT02-10
57	2mm extension pressure tube 15cm with 3 way stopcock red	CM-R3-PT02-15
58	2mm extension pressure tube 30cm with 3 way stopcock red	CM-R3-PT02-30
59	2mm extension pressure tube 140cm with 3 way stopcock red	CM-R3-PT02-140
60	3mm extension pressure tube 5cm with 3 way stopcock red	CM-R3-PT03-05
61	3mm extension pressure tube 10cm with 3 way stopcock red	CM-R3-PT03-10
62	3mm extension pressure tube 15cm with 3 way stopcock red	CM-R3-PT03-15

NO	NAME OF DEVICE	IDENTIFIER
63	3mm extension pressure tube 30cm with 3 way stopcock red	CM-R3-PT03-30
64	3mm extension pressure tube 140cm with 3 way stopcock red	CM-R3-PT03-140
65	3mm extension pressure tube (red line) 5cm with 3 way stopcock red	CM-R3-PT03-05R
66	3mm extension pressure tube (red line) 10cm with 3 way stopcock red	CM-R3-PT03-10R
67	3mm extension pressure tube (red line) 15cm with 3 way stopcock red	CM-R3-PT03-15R
68	3mm extension pressure tube (red line) 30cm with 3 way stopcock red	CM-R3-PT03-30R
69	3mm extension pressure tube (red line) 140cm with 3 way stopcock red	CM-R3-PT03-140R
70	4mm extension pressure tube 5cm with 3 way stopcock red	CM-R3-PT04-05
71	4mm extension pressure tube 10cm with 3 way stopcock red	CM-R3-PT04-10
72	4mm extension pressure tube 15cm with 3 way stopcock red	CM-R3-PT04-15
73	4mm extension pressure tube 30cm with 3 way stopcock red	CM-R3-PT04-30
74	4mm extension pressure tube 140cm with 3 way stopcock red	CM-R3-PT04-140
75	4mm extension pressure tube (red line) 5cm with 3 way stopcock red	CM-R3-PT04-05R
76	4mm extension pressure tube (red line) 10cm with 3 way stopcock red	CM-R3-PT04-10R
77	4mm extension pressure tube (red line) 15cm with 3 way stopcock red	CM-R3-PT04-15R
78	4mm extension pressure tube (red line) 30cm with 3 way stopcock red	CM-R3-PT04-30R
79	4mm extension pressure tube (red line) 140cm with 3 way stopcock red	CM-R3-PT04-140R
80	4mm Extension Tube 5cm with 3 way stopcock yellow	CM-Y3-ET04-05
81	4mm Extension Tube 10cm with 3 way stopcock yellow	CM-Y3-ET04-10
82	4mm Extension Tube 15cm with 3 way stopcock yellow	CM-Y3-ET04-15
83	4mm Extension Tube 30cm with 3 way stopcock yellow	CM-Y3-ET04-30
84	4mm Extension Tube 140cm with 3 way stopcock yellow	CM-Y3-ET04-140
85	4mm Extension Tube (yellow line) 5cm with 3 way stopcock yellow	CM-Y3-ET04-05Y
86	4mm Extension Tube (yellow line) 10cm with 3 way stopcock yellow	CM-Y3-ET04-10Y
87	4mm Extension Tube (yellow line) 15cm with 3 way stopcock yellow	CM-Y3-ET04-15Y
88	4mm Extension Tube (yellow line) 30cm with 3 way stopcock yellow	CM-Y3-ET04-30Y

NO	NAME OF DEVICE	IDENTIFIER
89	4mm Extension Tube (yellow line) 140cm with 3 way stopcock yellow	CM-Y3-ET04-140Y
90	4mm Extension Tube 5cm for Infusion Therapy with 3way stopcock blue	CM-3WB-ET04-05
91	4mm Extension Tube 10cm for Infusion Therapy with 3way stopcock blue	CM-3WB-ET04-10
92	4mm Extension Tube 15cm for Infusion Therapy with 3way stopcock blue	CM-3WB-ET04-15
93	4mm Extension Tube 30cm for Infusion Therapy with 3way stopcock blue	CM-3WB-ET04-30
94	4mm Extension Tube 140cm for Infusion Therapy with 3way stopcock blue	CM-3WB-ET04-140
95	4mm Extension Tube (blue line) 5cm for Infusion Therapy with 3way stopcock blue	CM-3WB-ET04-05B
96	4mm Extension Tube (blue line) 10cm for Infusion Therapy with 3way stopcock blue	CM-3WB-ET04-10B
97	4mm Extension Tube (blue line) 15cm for Infusion Therapy with 3way stopcock blue	CM-3WB-ET04-15B
98	4mm Extension Tube (blue line) 30cm for Infusion Therapy with 3way stopcock blue	CM-3WB-ET04-30B
99	4mm Extension Tube (blue line) 140cm for Infusion Therapy with 3way stopcock blue	CM-3WB-ET04-140B
100	3mm extension pressure tube 5cm with 3 way stopcock red	CM-3WR-PT03-05
101	3mm extension pressure tube 10cm with 3 way stopcock red	CM-3WR-PT03-10
102	3mm extension pressure tube 15cm with 3 way stopcock red	CM-3WR-PT03-15
103	3mm extension pressure tube 30cm with 3 way stopcock red	CM-3WR-PT03-30
104	3mm extension pressure tube 140cm with 3 way stopcock red	CM-3WR-PT03-140
105	4mm Extension Tube 5cm for Infusion Therapy with 3way stopcock red	CM-3WR-PT04-05
106	4mm Extension Tube 10cm for Infusion Therapy with 3way stopcock red	CM-3WR-PT04-10
107	4mm Extension Tube 15cm for Infusion Therapy with 3way stopcock red	CM-3WR-PT04-15
108	4mm Extension Tube 30cm for Infusion Therapy with 3way stopcock red	CM-3WR-PT04-30
109	4mm Extension Tube 140cm for Infusion Therapy with 3way stopcock red	CM-3WR-PT04-140
110	3mm extension pressure tube (red line) 5cm with 3 way stopcock red	CM-3WR-PT03-05R
111	3mm extension pressure tube (red line) 10cm with 3 way stopcock red	CM-3WR-PT03-10R

NO	NAME OF DEVICE	IDENTIFIER
112	3mm extension pressure tube (red line) 15cm with 3 way stopcock red	CM-3WR-PT03-15R
113	3mm extension pressure tube (red line) 30cm with 3 way stopcock red	CM-3WR-PT03-30R
114	3mm extension pressure tube (red line) 140cm with 3 way stopcock red	CM-3WR-PT03-140R
115	4mm Extension Tube (red line) 5cm for Infusion Therapy with 3way stopcock red	CM-3WR-PT04-05R
116	4mm Extension Tube (red line) 10cm for Infusion Therapy with 3way stopcock red	CM-3WR-PT04-10R
117	4mm Extension Tube (red line) 15cm for Infusion Therapy with 3way stopcock red	CM-3WR-PT04-15R
118	4mm Extension Tube (red line) 30cm for Infusion Therapy with 3way stopcock red	CM-3WR-PT04-30R
119	4mm Extension Tube (red line) 140cm for Infusion Therapy with 3way stopcock red	CM-3WR-PT04-140R
120	4mm Extension Tube 5cm for Infusion Therapy with 3way stopcock white	CM-3WW-ET04-05
121	4mm Extension Tube 10cm for Infusion Therapy with 3way stopcock white	CM-3WW-ET04-10
122	4mm Extension Tube 15cm for Infusion Therapy with 3way stopcock white	CM-3WW-ET04-15
123	4mm Extension Tube 30cm for Infusion Therapy with 3way stopcock white	CM-3WW-ET04-30
124	4mm Extension Tube 140cm for Infusion Therapy with 3way stopcock white	CM-3WW-ET04-140
125	4mm Extension Tube (white line) 5cm for Infusion Therapy with 3way stopcock white	CM-3WW-ET04-05W
126	4mm Extension Tube (white line) 10cm for Infusion Therapy with 3way stopcock white	CM-3WW-ET04-10W
127	4mm Extension Tube (white line) 15cm for Infusion Therapy with 3way stopcock white	CM-3WW-ET04-15W
128	4mm Extension Tube (white line) 30cm for Infusion Therapy with 3way stopcock white	CM-3WW-ET04-30W
129	4mm Extension Tube (white line) 140cm for Infusion Therapy with 3way stopcock white	CM-3WW-ET04-140W
130	4mm Extension Tube 5cm for Infusion Therapy with 3way stopcock yellow	CM-3WY-ET04-05

NO	NAME OF DEVICE	IDENTIFIER
131	4mm Extension Tube 10cm for Infusion Therapy with 3way stopcock yellow	CM-3WY-ET04-10
132	4mm Extension Tube 15cm for Infusion Therapy with 3way stopcock yellow	CM-3WY-ET04-15
133	4mm Extension Tube 30cm for Infusion Therapy with 3way stopcock yellow	CM-3WY-ET04-30
134	4mm Extension Tube 140cm for Infusion Therapy with 3way stopcock yellow	CM-3WY-ET04-140
135	4mm Extension Tube (yellow line) 5cm for Infusion Therapy with 3way stopcock yellow	CM-3WY-ET04-05Y
136	4mm Extension Tube (yellow line) 10cm for Infusion Therapy with 3way stopcock yellow	CM-3WY-ET04-10Y
137	4mm Extension Tube (yellow line) 15cm for Infusion Therapy with 3way stopcock yellow	CM-3WY-ET04-15Y
138	4mm Extension Tube (yellow line) 30cm for Infusion Therapy with 3way stopcock yellow	CM-3WY-ET04-30Y
139	4mm Extension Tube (yellow line) 140cm for Infusion Therapy with 3way stopcock yellow	CM-3WY-ET04-140Y
140	4mm Extension Tube 5cm with 3 way stopcock Needle Free Valve	CM-NF-B3-ET04-05
141	4mm Extension Tube 10cm with 3 way stopcock Needle Free Valve	CM-NF-B3-ET04-10
142	4mm Extension Tube 15cm with 3 way stopcock Needle Free Valve	CM-NF-B3-ET04-15
143	4mm Extension Tube 30cm with 3 way stopcock Needle Free Valve	CM-NF-B3-ET04-30
144	4mm Extension Tube 140cm with 3 way stopcock Needle Free Valve	CM-NF-B3-ET04-140