



DEVICE DETAILS

NAME OF DEVICE	Aneroid Sphygmomanometer (Aneroid Gauge Sphygmomanometer)
ESTABLISHMENT NAME	CYPRESS MEDIC SDN BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GA8626619-31311
BRAND NAME	Rossmax
MEDICAL DEVICE CATEGORY	HEALTHCARE FACILITY PRODUCTS AND ADAPTATIONS
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	Aneroid Sphygmomanometer is made of Gauge, Bulb, Bladder, Cuff, Stethoscope earpiece, Stethoscope pipe and Stethoscope tip etc. Sphygmomanometer adopts European indicating tray, the minimum indication value is 2mm Hg; the hand rotates calmly, it does not jump or do not move. The outer case is made from AL.alloy or Zn. Alloy and backed by using all kinds of thin pattern lacquer. When cuff is wrapped, it does not cause stimulation, allergic reaction and other harm.
DEVICE INTENDED PURPOSE	It is used to measure the blood pressure of human body.
VALIDITY DATE OF REGISTRATION	22/07/2019 - 21/07/2024

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Aneroid Sphygmomanometer (Aneroid Gauge Sphygmomanometer)	GD101
2	Aneroid Sphygmomanometer (Aneroid Gauge Sphygmomanometer)	GD102
3	Aneroid Sphygmomanometer (Aneroid Gauge Sphygmomanometer)	GD111
4	Aneroid Sphygmomanometer (Aneroid Gauge Sphygmomanometer)	GD112
5	Aneroid Sphygmomanometer (Aneroid Gauge Sphygmomanometer)	GB101
6	Aneroid Sphygmomanometer (Aneroid Gauge Sphygmomanometer)	GB102
7	Aneroid Sphygmomanometer (Aneroid Gauge Sphygmomanometer)	GB111
8	Aneroid Sphygmomanometer (Aneroid Gauge Sphygmomanometer)	GB112