



DEVICE DETAILS

NAME OF DEVICE	MEDISANA PULSE OXIMETER PM 100 CONNECT
ESTABLISHMENT NAME	ZUELLIG PHARMA SDN.BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GC9976721-74944
BRAND NAME	MEDISANA
MEDICAL DEVICE CATEGORY	MD 1100 - GENERAL ACTIVE MEDICAL DEVICES
DEVICE GROUPING TYPE	SINGLE
DEVICE DESCRIPTION	<p>The Pulse Oximeter PM 100 CONNECT is a portable non-invasive device used to measure and record the oxygen saturation of arterial hemoglobin (SpO₂) and pulse rate of users at any moment. It is a small and compact device which is applicable in the hospital and home use with simple and convenient operation. However, it is reminded that the Pulse Oximeter PM 100 CONNECT is not suitable for continuous blood oxygen monitoring. The measuring range for the SpO₂ is 70% to 99% with 1% resolution while the measuring range for the pulse rate is 30 to 235 beats per minute (bpm) with a resolution of 1 bpm.</p>
DEVICE INTENDED PURPOSE	<p>The Pulse Oximeter PM 100 CONNECT is intended for spot-checking of oxygen saturation of arterial hemoglobin (SpO₂) and pulse rate users. It is not suitable for continuous monitoring.</p>
VALIDITY DATE OF REGISTRATION	20/09/2021 - 19/09/2026

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	MEDISANA PULSE OXIMETER PM 100 CONNECT	48506